Restek Electronic Leak Detector

(Cat.# 28500)

Available Languages

Instruction Manual	English	Page 3
Gebrauchsanweisung	Deutsch	Seite 11
Notice d'utilisation	Français	Page 19
Manuale di istruzioni	Italiano	Pagina 27
Manual de instrucciones	Español	Página 35
使用说明书	Chinese	页 43
取扱説明書	Japanese	ページ 51

Table of Contents

- 1.0 Introduction 4
- 2.0 Battery Charging 5
- 3.0 Powering Up 5
- 4.0 Zeroing the Unit 5
- 5.0 Prior to Operation 6
- 6.0 Detecting Leaks 6
- 7.0 Specifications..... 7
- 8.0 Maintenance 7
- 9.0 Troubleshooting 7
- 10.0 Interpretation of Results 8
- 11.0 Back Label Explanation 9
- 12.0 Case Specifications 10
- 13.0 Service 10
- Declaration of Conformity 59

Restek Electronic Leak Detector

Operating Instructions

1.0 Introduction

Restek's portable electronic leak detector is specifically designed for use with gas chromatography (GC) systems. It detects small leaks of any gas with a thermal conductivity different from air. The reference gas inlet (Figure 1) draws in ambient air for comparison to gas drawn into the sample probe. The intensity of a leak is indicated by both an LED light display and an audible alarm. The following table explains how to interpret the LED display.

NOTE: The leak detector is not a quantitative device, it is a qualitative measure of the location, intensity, and nature of a gas leak.

It is best practice to use a leak detector daily to check critical seals (for example, septa, column nuts, reducing nuts, and fittings on gas lines). In addition to being used regularly, the leak detector should be charged regularly to ensure the battery pack continues to hold its charge.

If this instrument is used in any manner other than described in the manual, the CE declaration is void.

The leak detector is manufactured by Restek and is backed by our 100% satisfaction guarantee.

Should you require assistance at anytime regarding our leak detector, please contact Restek Customer Service.

Leak Detector LED Indicator Display and Interpretation				
	Unit State - Powered On/Off	LED Color	LED State	Interpretation
Unit Battery State	On	Blue > Red & Yellow	On/Blue flashes once, red and yellow cycle for 15 seconds/Off	Unit is powering on and running startup sequence
	On	Blue	On/Steady	Ready for use, battery does not require charging
	On	Blue	On/Constant Flashing	IMPORTANT: Unit must be charged
	Power Button Depressed	NONE	NONE	Battery has fully discharged, charge immediately
Unit Charging State	On/Off	Green	On/Flashing	Unit is plugged in to the AC adaptor or USB cable and is charging a very low battery (trickle charging)
	On/Off	Green	On/Steady	Unit is plugged in to AC adaptor or USB cable and is in full charge mode
	On/Off	Green	Off	Unit is plugged into AC adaptor or USB cable and the unit is fully charged
Unit Operation	On	Red or Yellow	On/Steady	Unit is indicating a difference in thermal conductivity between probe and reference
	On	Red & Yellow	On/Cycling for 4 seconds/Off	Zero set button has been depressed and unit is rezeroing

2.0 Battery Charging

The leak detector should be fully charged prior to use. Only use the AC adaptor or USB cable, both of which are provided with the unit. To charge the battery using the AC adaptor, first install the correct plug for your country's AC outlets onto the provided AC adaptor. Insert the AC adaptor into an electrical outlet, and then insert the barrel plug on the other end of the AC adaptor into the connector on the bottom of the leak detector unit. To charge the battery using the USB cable, simply plug the USB cable into the unit and into a USB port. While charging, the green battery charge indicator LED will illuminate. When the battery is fully charged, the green battery charge indicator LED will no longer be illuminated.

When the leak detector's charge is low, the blue LED located between the red and yellow LEDs will begin to flash. It is important to charge the leak detector when the blue light is flashing to avoid complete battery discharge and possible damage to the battery.

NOTE: If the battery is fully discharged, no LED will illuminate.

CAUTION: For the health of the battery, it is not recommended to continuously operate the leak detector when it is plugged in and charging. However, if your battery has been allowed to fully discharge and will no longer hold a charge, the unit can be operated while plugged in receiving charge directly from the power supply. To operate the leak detector while it is plugged in, follow the instructions for regular operation, and when the blue LED illuminates steady or flashing the unit is ready to use.

CAUTION: The leak detector contains a nickel-metal hydride (NiMH) battery. Like other battery-operated devices, if the battery is left for long periods of time without being charged, it can discharge to a point where the battery cannot accept a charge. **We recommend that you charge the battery at least once every 3 months. If your leak detector is used often for routine maintenance, as recommended, more frequent charging may be necessary.**

CAUTION: Replacement of the rechargeable battery in this unit is performed at the factory. There are no serviceable parts inside this unit. Opening the case or tampering with the internal parts will void the factory warranty.

NOTE: Recharging a low battery can take approximately 6 hours.

NOTE: The leak detector may get warm to the touch during charging, but this is normal and not a safety concern.

3.0 Powering Up

Remove the probe assembly from the holder. Depress and hold the power button (Figure 1) until the unit responds by initiating a starting sequence that lasts approximately 15 seconds, during which the red and yellow LEDs will flash. When the red and yellow LEDs stop flashing and the blue LED illuminates, the unit is ready to use or zero. During the **starting sequence, DO NOT** attempt to zero the unit.

4.0 Zeroing the Unit

The instrument may need to be zeroed periodically between uses, especially if it is moved from room to room or between areas of differing temperature or humidity. **Do not attempt to zero the unit while the probe is stored in the holder.** The probe **MUST** be removed from the probe holder (Figure 1) before zeroing the unit. To re-zero, press the zero set button. The unit will run a re-zeroing sequence for approximately 4 seconds, during which the red and yellow LEDs will flash. When all LED lights stop flashing and the blue LED light is lit, the unit is ready for use.

NOTE: To avoid false readings, do not attempt to use or zero the unit while the zeroing sequence is in progress.


5.0 Prior to Operation

Verify the operation of the leak detector before each use by sampling gas from a GC split vent or other source of gas (other than air). Also, visually inspect the probe tip, reference gas inlet, and exhaust port for obstructions (Figure 1).

IMPORTANT: Fittings being checked must be clean and dry; liquids, dust, and other debris may damage the leak detector if drawn into the probe. Testing for leaks around hot fittings could increase the risk of burns. Hot fittings may also affect the operation of the leak detector.

NOTE: The leak detector responds to almost any gas. Solvent vapors, split vent exhaust, or even strong air currents around the probe or reference inlet can cause instability or false positive readings. Be careful not to breathe into the reference inlet when checking for leaks, and do not cover/block the reference inlet with your hand.

Figure 1: Leak detector schematic.

6.0 Detecting Leaks

With the leak detector turned on and zeroed, slowly move the probe tip around fittings and other potential leak sources. If the leak detector senses a gas other than air, the LEDs will begin to light, with more lights indicating a more significant leak. On the third red, or second yellow LED, an audible tone will begin to beep. The more LEDs that illuminate, the faster the beep frequency. When the last red or yellow LED illuminates, the beep becomes a steady tone. The red LEDs indicate a helium or hydrogen leak. The yellow LEDs indicate a nitrogen, argon, or carbon dioxide leak.

Remove the probe from the vicinity of the leak and allow the unit to return to zero. If a large amount of gas has entered the probe, it may take a few seconds for the instrument to clear that gas and reset itself. **Do not attempt to zero the unit while it is exhausting the gas from the probe. This may cause the unit to malfunction.** Once reset, place the probe near the leak again to confirm its location. **The reference gas inlet (Figure 1) must not be restricted or the unit will not operate correctly. Similarly, the exhaust port allows the gas being tested to exit the leak detector and must remain unobstructed.** The exhaust port is located in the probe holder.

CAUTION: This unit should only be used to detect trace amounts of hydrogen arising from a small leak in a non-flammable environment, e.g., laboratory room air, etc.

NOTE: To disable the audible beep during leak detection, depress and hold the zero set button for 2–3 seconds. After you hear a steady tone for 1 second, release the button; the beep function is disabled. To turn the beep function on again, depress and hold the zero set button. The beep function is always enabled at power-up.

NOTE: The leak detector will power down after 6 minutes of operation. This feature prevents excess battery discharge if the unit is accidentally left on.

7.0 Specifications

- Altitude: Up to 2000 m above mean sea level
- Ambient Temperature: 50–98.6 °F (10–37 °C)
- Battery Rating: 12 hours operation
- Certifications: CE (EU, Korea, Japan, Australia); UKCA; CSA/UL tested, not listed.
- Compliance: WEEE, CEC, China RoHS 2
- Humidity Range: 0–97%
- Indoor Use Only
- Manufacturer: Restek Corporation, 110 Benner Circle, Bellefonte, PA 16823
- Pollution Degree: 2
- Power Rating: 5 volts DC, 1 A (AC adaptor and USB charging cable supplied)
- Warranty: 1-year warranty



8.0 Maintenance

Avoid spilling liquids onto the unit or it may malfunction. If a liquid is spilled onto the unit, turn off the power immediately, remove excess liquids with a dry towel, and let the unit sit until the liquid dries. Dust and debris can enter the probe tip of the leak detector and, over time, can clog the small-bore tubing inside the unit. To prevent this, clean the probe tip periodically. To clean the probe tip, unscrew the cap to expose the brush (Figures 2 and 3). Gently clean the probe tip brush using a small cleaning brush or your fingers to remove dust and debris, then replace the cap. Do not use liquids to clean the probe. Liquids can damage the leak detector if drawn in through the probe.

Information on where to have the unit sent for maintenance or service* is listed at the end of this document.

9.0 Troubleshooting

Problem	Possible Cause	Suggested Solution
Change in sensitivity is observed	Probe clogged	Clean the probe tip to remove any debris
	Probe line punctured, damaged or plugged	Visually inspect probe line for holes*
Response decreased	Detector not zeroed	Re-zero detector
LEDs stays lit during operation	Detector re-zeroed before unit was completely reset	Allow adequate time for detector to vent gas, then re-zero
	Reference gas or probe inlet covered by hand or other object	Remove obstruction
Does not power up	Batteries need to be charged	Charge unit
Does not power up after charging	Batteries may be damaged	Unit can operate while plugged into adaptor or USB cable, but unit should be returned to Restek to have the battery replaced*

*Contact returns@restek.com for return instructions for servicing a damaged unit.


10.0 Interpretation of Results

Figure 4 illustrates the leak detector's LED light response range. A greater number of lit red or yellow LED lights correlates in general to a larger leak.

NOTE: The leak detector is not a quantitative device, rather it is designed to detect leaks in gas line connections commonly associated with laboratory equipment.

Gas	Minimum Detectable Leak Rate (atm cc / sec)	Indicating LED Light Color
Helium	1.0×10^{-5}	Red
Hydrogen**	1.0×10^{-5}	Red
Nitrogen	1.4×10^{-3}	Yellow
Argon	1.0×10^{-4}	Yellow
Carbon dioxide	1.0×10^{-4}	Yellow

Figure 4: LED light response chart for the leak detector. A 1:1 ratio (probe:reference) indicates no leak present. Red LED lights indicate the presence of helium and/or hydrogen. Yellow LED lights indicate the presence of nitrogen, argon, and/or carbon dioxide.

****CAUTION:** This unit should only be used to detect trace amounts of hydrogen arising from a small leak in a nonflammable environment, e.g., laboratory room air, etc.


Tip drift

Tip drift is the phenomenon when a false LED light response is registered as the probe or the unit is quickly moved. Tip drift is inherent to all dual thermistor leak detector technology and is based in large part on the asymmetry of the flow cells; shaking or tipping the unit influences the air flow profiles, which impacts the rates of heat exchange. If the device is functioning normally, the LED light signal will return to zero in 3–5 seconds after the unit is held still. In extreme cases, the unit may need to be re-zeroed before using. To avoid tip drift, be sure to hold the unit body steady while making measurements with the probe tip.


11.0 Back Label Explanation

- 1. Product name.
- 2. Product catalog number.
- 3. Product serial number.
- 4. Passed California Energy Commission (CEC) Battery Charging System (BCS) testing.
- 5. CE mark: see Declaration of Conformity
- 6. KC mark - Korea Certification
- 7. Electrical parameters.
- 8. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
 - 1) This device may not cause harmful interference.
 - 2) This device must accept any interference received, including interference that may cause undesired operation.
- 9. Manufacturer company name, address, and contact information.
- 10. AC adaptor polarity
- 11. VCCI mark - VCCI certification
- 12. RCM (Regulatory Compliance Mark) - Australia certification
- 13. China RoHS 2
- 14. Unit is WEEE compliant.
- 15. Unit is UKCA compliant.
- 16. This Class A digital apparatus complies with Canadian ICES-003.
- 17. Units must be sent back to Restek Corporation for service.

For the most up-to-date information, see our website www.restek.com/leakdetector

12.0 Case Specifications

13.0 Service

The Restek leak detector carries a one-year limited warranty from time of purchase. Please have the leak detector serial number available when calling Restek with any concerns you may have.

For questions, problems, or repair services:

Within the U.S.:

Call Restek Customer Service at
1-800-356-1688 or 1-814-353-1300, ext. 3.

Outside the U.S.:

Contact your local
Restek representative.

Optional Accessories

Soft-Side Carry/Storage Case

Ideal for storing the leak detector in smaller spaces such as a toolbox.
cat.# 22657

Small Probe Adaptor

Verify hard-to-reach leaks using the small probe adaptor.
cat.# 22658

Charging Kit

Includes a universal power adaptor and USB charging cable.
cat.# 28502

USB Charging Cable

For convenient USB port charging.
cat.# 28501

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DECLARATION OF CONFORMITY

Application of Council Directive (s):	EMC Directive 2014/30/EU Safety Low Voltage Directive 2014/35/EU RoHS Directive 2011/65/EU
Standards to Which Conformity is Declared:	EN 61326-1:2013 FCC CFR 47, Part 15, subpart B:2017, Class A ICES-003, Issue 6:2016, Class A ITE VCCI/V-3:2015.04 IEC 61326-1 Ed. 2 (2012) / AS/NZS CISPR 11:2017 KN 32:2015 / KN 35:2015 IEC 61010-1:2010, AMD1:2016- Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements IEC 61010-1:2010 CANADA / US NATIONAL DIFFERENCES National standard CSA C22.2 No. 61010-1-12 / UL 61010-1:2012 IEC 61010-1:2010 SWITZERLAND NATIONAL DIFFERENCES SN EN 61010-1:2010 IEC 61010-1:2010 JAPAN NATIONAL DIFFERENCES IEC 61010-1:2010 Ed. 3
Declarer's Name:	Restek Corporation
Declarer's Address:	110 Benner Circle Bellefonte, PA 16823
Type of Equipment:	Leak Detector 6
Catalog Number(s):	28500
Person responsible for making this Declaration:	Gerald Johnston
Position Title:	Research & Development Engineer III
Date:	05/22/2019
Signature:	

Questions about this or any other Restek product? Contact us or your local Restek representative (www.restek.com/contact-us).

Restek patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2022 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com

#204-07-014 Rev. date: 07/22